


Instruction Manual

Elbow Crutch




FCTAL02

SectorCare

<https://sectorcare.com.au>



Product Installation



Parts List



Upper Shaft*1 Lower Shaft*1

1. Take out all parts and accessories after opening the package.
2. Align the connecting joints of the upper Shaft and Lower Shaft, insert into place, then secure and tighten.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.