

Instruction Manual

Toilet Aid with Leg



TSCAL01

SectorCare

<https://sectorcare.com.au>



Product Installation



01 Armrest

02 Toilet Seat

03 Adjustable knob

04 Height Adjustable Leg

05 Anti-Slip Foot Cap

Parts List:



Commode Seat*1



Metal Leg Tube*4



Armrest Grip*2



U-shaped Clip*4


1. Take out all products and accessories after opening the package.



2. Place the commode seat upside down. Align the four leg tubes with the four connectors respectively, insert them and secure with U-shaped clips, as shown in the figure.



3. Align the armrest with the connector on the leg tube, insert and secure it, as shown in the figure.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

© 2026 FRG Pty Ltd. All rights reserved.